SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The SSCP is not intended to replace the Instructions For Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

The following information is intended for users/healthcare professionals.

Following this information there is a summary intended for patients.

Manufacturer's reference number for the SSCP - 8032472BIOANAbY2AA_SSCP

1. DEVICE IDENTIFICATION AND GENERAL INFORMATION

1.1 Device trade name(s)

NEOBIANACID, NEOBIANACID PEDIATRIC

1.2 Manufacturer's name and address

Aboca S.p.A. Società Agricola - Loc. Aboca 20 – 52037 Sansepolcro (AR) – Italy

1.3 Manufacturer's single registration number (SRN)

IT-MF-000029499

1.4 Basic UDI-DI

8032472BIOANAbY2AA

1.5 Medical device nomenclature description / text

G0401 - ORALLY ADMINISTERED DEVICES FOR THE THERAPY OF GASTRO-INTESTINAL DISORDERS

1.6 Class of device

Class III

1.7 Year when the first certificate (CE) was issued covering the device

2014 (certification according to Directive 93/42/EEC)

1.8 Authorised representative if applicable; name and the SRN

Not applicable

1.9 NB's name (the NB that will validate the SSCP) and the NB's single identification number

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2. INTENDED USE OF THE DEVICE

2.1 Intended purpose

Neobianacid is indicated for the treatment of gastroesophageal reflux and stomach pain (painful dyspepsia), characterized by acid regurgitation, burning, pain, irritative cough, dysphonia.

Neobianacid is also indicated in the treatment of symptoms associated with digestive difficulties (non-painful dyspepsia), such as acidity, heaviness, and bloating after meals, both recurrent and in cases of copious meals.

Neobianacid can be useful in prevention, in situations that cause irritation of the mucosa (e.g., taking anti-inflammatory drugs, periods of mental and physical stress).

The product is also suitable for use in pregnancy, during lactation and in children over 6 years of age.

Neobianacid Pediatric is indicated from one year of age for the treatment of gastroesophageal reflux and acidity characterized by regurgitation, vomiting, irritative cough and abdominal pain.

Neobianacid Pediatric is also indicated in the treatment of symptoms associated with dyspepsia (digestive difficulty), such as pain and heartburn, heaviness and bloating after meals.

2.2 Indication(s) and target population(s)

Neobianacid is indicated in adults and children over 6 years of age. The device is also suitable for use in pregnancy and during breastfeeding.

Neobianacid Pediatric is indicated in children over 1 year of age.

2.3 Contraindications and/or limitations

The device shall not be used in case of individual hypersensitivity or allergy to one or more components.

However, if symptoms persist, consult your doctor or pharmacist.

It is recommended to take the device at least two hours before or after taking medicines, particularly tetracyclines. There are no known interactions of Neobianacid or Neobianacid Pediatric with other medical devices or foods.

3. DEVICE DESCRIPTION

3.1 Description of the device

Neobianacid is a medical device available in tablets and granular sachets. Neobianacid Pediatric is available in granular sachets.

Neobianacid needs to be dissolved in the mouth and Neobianacid Pediatric needs to be dissolved in little water and then ingested.

Neobianacid and Neobianacid pediatric act rapidly by forming a barrier effect film with the following activities:

PROTECTIVE

Thanks to its high affinity for the gastroesophageal mucosa, it forms a film that protects and promotes the physiological reparative processes.

LOCAL NEUTRALIZER OF ACIDITY

Strengthens the local neutralizing properties of gastric mucus, hindering the contact of acids with the gastric mucosa and preventing its damage.

ANTIOXIDANT

Counteracts irritation due to free radicals. This way it protects the mucosa from the contact with gastric juices and irritating substances.

Neobianacid and Neobianacid pediatric modulate acidity without suppressing it, thus respecting the physiology of the stomach and the functionality of the intestinal flora.

Neobianacid and Neobianacid pediatric are 100% natural, biodegradable and with organic ingredients. They contain Poliprotect, a vegetal and mineral complex. They do not contain synthetic, semi-synthetic or genetically modified substances and the manufacturing processes with which it is obtained does not use synthetic solvents or other chemically synthesized substances.

The qualitative composition of Neobianacid is described below:

One 1.55 g tablet or sachet contains:

- Poliprotect® natural complex 528 mg, based on
 - Polysaccharides (LPME extractive fractions of Mallow* and Marshmallow*; dehydrated Aloe vera leaf gel * - 105 mg) and
 - Minerals (Limestone and Nahcolite 423 mg),
- LPME extractive flavonoid fraction of Chamomile* and Liquorice*, - 46 mg;

Brown Sugar*, Arabic Gum, Natural mint flavour or Natural lemon flavour.

*Ingredient from organic farming.

Each tablet contains 136 mg (3.4 mmol) of calcium.

The qualitative composition of Neobianacid pediatric is described below:

A 775 mg sachet of Neobianacid pediatric contains:

- Poliprotect® natural complex 264 mg, based on
 - Polysaccharides (LPME extractive fractions of Mallow* and Marshmallow *; dehydrated Aloe vera leaf gel* - 52 mg) and
 - Minerals (Limestone and Nahcolite 212 mg);
- LPME extractive flavonoid fraction of Chamomile* and Liquorice* - 23 mg;

Brown sugar*; Arabic gum; Natural flavour.

*Ingredient from organic farming.

Each sachet contains 68 mg (1.7 mmol) of calcium.

In both Neobianacid and Neobianacid pediatric, the extractive fractions contained in the product are obtained through an exclusive technological process, called LPME (Liquid-Phase Microextraction), developed by Aboca research.

3.2 A reference to previous generation(s) or variants if such exist, and a description of the differences

Neobianacid has been certified according to Directive 93/42/EEC since 2014.

The device is available in different packaging sizes in the following dosage forms:

- Tablets to be dissolved in the mouth (adults and children starting from 6 years of age)
- Granular sachets to be dissolved in the mouth (adults and children starting from 6 years of age)
- Pediatric variant in granular sachets to be dissolved in a little water (children starting from 1 year of age)

Neobianacid in the dosage forms for adults and children starting from 6 years of age can be available in natural mint flavour or natural lemon flavour.

Neobianacid pediatric is available in strawberry taste.

3.3 Description of any accessories which are intended to be used in combination with the device

Not applicable – no accessories are intended to be used in combination with Neobianacid

3.4 Description of any other devices and products which are intended to be used in combination with the device

Not applicable as there are no devices and products intended to be used in combination with Neobianacid. Nevertheless, Neobianacid can be used in combination with antisecretory drugs, specifically proton pump inhibitors (PPIs) and H2 antagonists, as well as during periods of suspension of treatment, as advised by a doctor.

4. RISKS AND WARNINGS

4.1 Residual risks and undesirable effects

The main residual risks identified for Neobianacid are the following:

- Subjective allergic reaction of the user – although the product is non-sensitising, it is possible that a small number of people may be allergic to one or more components of the product; the labelling and the instructions for use list the components of the product, so the user has all the information needed in order to verify the ingredients and therefore the safe use in cases of allergies.
- Microbial contamination in case of poor storage and/or use of the product – in order to minimize this risk, the storage conditions of the product are clearly reported on the labelling and instructions for use.

Neobianacid is normally well tolerated, however, if the user notices an adverse effect potentially attributable to the product or its use, it is recommended to inform a doctor and/or pharmacist and the Manufacturer Aboca.

4.2 Warnings and precautions

Do not use in case of individual hypersensitivity or allergy to one or more components. If symptoms persist, seek medical advice. The glycyrrhizin content from the Liquorice in Neobianacid is lower than that reported to cause a possible increase in blood pressure.

4.3 Other relevant aspects of safety, including a summary of any field safety corrective action (FSCA including FSN) if applicable

No FSCA (Field Safety Corrective Action) has been implemented for the device Neobianacid since its first date of commercialization in 2014.

5. SUMMARY OF CLINICAL EVALUATION AND POST-MARKET CLINICAL FOLLOW-UP (PMCF)

The efficacy and safety of Neobianacid are demonstrated by clinical data coming from:

- two randomized clinical trials;
- several large-scale surveys conducted in a real-world setting addressed to both healthcare professionals and patients;
- safety data collected from the manufacturer's Vigilance system.

5.1 Summary of clinical data related to equivalent device, if applicable

Not applicable to Neobianacid

5.2 Summary of clinical data from conducted investigations of the device before the CE-marking, if applicable

Not applicable since it is a legacy medical device

5.3 Summary of clinical data from other sources

As part of the post-market surveillance activities, two randomized, double-blind, parallel group clinical trials were conducted on Neobianacid. Also other data on the product gathered in a real world setting are available and are reported below.

5.3.1 Clinical trials

About the clinical trials, one is published, it is a non-inferiority trial conducted on Neobianacid to demonstrate its clinical efficacy and safety compared to Omeprazole 20 mg/day in the treatment of non-erosive gastroesophageal reflux disease and painful dyspepsia.

The other clinical trial aimed to evaluate the clinical efficacy of Neobianacid compared to its equivalent medical device Mipolixin (manufactured by Aboca) in improving the overall symptoms severity of functional dyspepsia and/or heartburn.

A summary of the two clinical trials is reported below.

- **Published Clinical trial on the American Journal of Gastroenterology: "Poliprotect vs Omeprazole in the Relief of Heartburn, Epigastric Pain, and Burning in Patients Without Erosive Esophagitis and Gastroduodenal Lesions: A Randomized, Controlled Trial"¹**

This was a randomized, double blind, double dummy, reference product controlled, parallel group, non-inferiority clinical study that compared the clinical efficacy of two different agents, the medical device

¹ Enrico Stefano Corazziari et al. Poliprotect vs Omeprazole in the Relief of Heartburn, Epigastric Pain, and Burning in Patients Without Erosive Esophagitis and Gastroduodenal Lesions: A Randomized, Controlled Trial. Am J Gastroenterol 2023;00:1–11. <https://doi.org/10.14309/ajg.0000000000002360>

Neobianacid and Omeprazole, in the relief of symptoms of heartburn and epigastric pain in subjects with non-erosive gastroesophageal reflux disease (NERD) and/or epigastric pain syndrome (EPS).

Study objectives

Primary objective of this study was to assess whether Neobianacid has a comparable symptomatic effect to Omeprazole on heartburn and epigastric pain in the treatment of adult NERD and EPS patients. Secondary objectives were to assess how fast Neobianacid provides symptomatic relief after the first assumption and to what extent might provide relief for other gastrointestinal symptoms. In addition, the study assessed the patient's degree of treatment satisfaction and the effects on quality of life, as well as the treatment's safety and tolerability. Exploratory aim of the study was to assess the gut microbiota change after Neobianacid and Omeprazole treatment.

Inclusion and exclusion criteria

Patients presenting heartburn diagnosed with NERD according to Montreal criteria and/or diagnosed with EPS according to Rome III criteria, were eligible to enter this study. Additional inclusion criteria were the following: age from 18 to 75 years old of both genders; absence of macroscopic lesions of gastro-esophageal mucosa at endoscopy performed during the screening/wash-out period if not already performed during the previous 3 years and the symptomatology was unchanged; a visual analogue scale (VAS) score for heartburn and/or epigastric pain between 30 and 70 on at least 6 out of 14 days of the screening/wash-out period; agreement not to change their alimentary habits and lifestyle throughout the study period. Main exclusion criteria were: history of no response of heartburn/epigastric pain to antacids or PPI; patients using any drug that could affect symptoms, such Non-steroidal anti-inflammatory drugs (NSAIDs), antacids, anti-secretory agents; antibiotic use within 2 months preceding the screening period; Barrett's oesophagus; history of gastrointestinal surgery; atopy; gallstones; persistent or recurring vomiting; anorexia; dysphagia or odynophagia; any chronic disease that may significantly affect the gastrointestinal upper tract. Pregnant or fertile women without contraception were also excluded. Moreover, patients were required to avoid any of the above-mentioned drug since the beginning of screening/wash-out period until the end of the study.

Study design and methods

The study design consisted of a 2-week screening/wash-out period, a double-blind 4-week randomized treatment period, followed by an open label 4-week period, with 5 visits in total. After the screening period, 275 patients fulfilling selection criteria were 1:1 randomised to the assigned treatment group for a double-blind 4-week treatment period. 141 patients (51.3% of randomised) received Neobianacid and 134 (48.7%) received Omeprazole. There were no important differences between groups in demographic characteristics. The study population included more females than males in both treatment groups.

The two study groups were submitted to 2 different regimens:

From Day 0 to Day 13

- GROUP A: Neobianacid five times daily + Omeprazole placebo once daily
- GROUP B: Omeprazole 20 mg once daily + Neobianacid placebo five times daily

From Day 14 to Day 27

- GROUP A: Neobianacid on demand + Omeprazole placebo once daily
- GROUP B: Omeprazole 20 mg once daily + Neobianacid placebo on demand

This blinded 4-week randomized treatment period was followed by the last open label 4 weeks of the study (from day 28 to day 55), during which all patients took Neobianacid only and on-demand, irrespective of previous treatment group.

A daily diary was to be completed at home for the whole duration of the study. An antacid (magaldrate oral gel) was allowed as rescue medication at any time during the study and recorded.

Endpoints

The primary efficacy endpoint of the study was the comparison between the 2 groups of heartburn and/or epigastric pain severity from baseline to day 13, as measured by means of a 100mm VAS (from “no symptoms” equal to 0mm to “overwhelming symptoms” equal to 100mm) reported on patients’ daily diary. Secondary efficacy endpoints were the following: comparison between the 2 groups of VAS score change for heartburn or epigastric pain severity at day 1, day 3 and day 7; comparison between groups of VAS score change for heartburn or epigastric pain severity from day 14 to day 27; difference between groups of quantity, number of days of use and proportion of patients using rescue medication during each phase of treatment period; between-group difference of on-demand intake of Neobianacid; changes in Gastrointestinal Quality of Life Index (GIQLI) and Gastrointestinal Symptom Rating Scale (GSRS) scores at day 14, day 28 and day 56 versus baseline between groups; Overall Treatment Evaluation (OTE) scores at day 14, day 28, day 56, between groups.

There was also an exploratory endpoint consisting in the characterization of the whole intestinal microbiota carried out by means of the analysis of faecal samples.

The safety and tolerability assessment included incidence of adverse events (AEs), treatment-related AEs and data from the laboratory test results.


Summary of results

The results of the primary efficacy endpoint relief of heartburn or epigastric pain severity (measured using a 100 mm VAS) assessed at Day 14 showed a comparable decrease in mean VAS score in both treatment groups (Figure 1).

The clinical benefit of Neobianacid and Omeprazole was evident as early as the first day of treatment and there were no significant differences between groups in the VAS score change 1, 3 and 7 days after treatment start.

The effects on heartburn or epigastric pain severity in the two groups assessed at Day 14 were maintained throughout the study, as assessed at Day 28 and Day 56, as shown in Figure 1.

Figure 1

Despite comparable on-demand intake of Neobianacid and VAS score changes between groups, from day 14 to day 27 the number of sachets of rescue medication used was significantly lower in Neobianacid group and remained significantly lower from day 28 to day 55 as well.

During the last unblinded 4 weeks of the study, the 2 groups did not differ in Neobianacid consumption and symptomatic benefit, which, in the omeprazole arm, was obtained alongside a significant increase in antacid rescue medicine as compared to the Neobianacid arm. This finding is in line with the notion that PPI suspension can be followed by symptoms worsening due to the gastric acid rebound effect. On note, rescue medication use in the Omeprazole group remained steady after shifting to on-demand Neobianacid intake (i.e, after daily Omeprazole intake was stopped) indicating that Neobianacid on-demand counteracted the predictable worsening of symptoms that follows the suspension of PPI treatment.

Consistently with the results of heartburn or epigastric pain severity, improvements from baseline to day 14, day 28 and day 56 in symptomatology and overall quality of life were observed in both treatment groups. Abdominal pain, diarrhoea and indigestion score changes as measured by means of Gastrointestinal Symptom Rating Scale (GSRS) did not significantly differ between the two groups at all time points.

The results indicate that Neobianacid, unlike Omeprazole, does not produce a breakdown of the gastric acid antimicrobial barrier. Thus, the study demonstrated a lower impact of Neobianacid on the intestinal microbiota, preserving the natural acid barrier of the stomach.

Overall, Neobianacid was very well tolerated. The proportion of patients with adverse events (AEs) was comparable in the two treatment groups and none of patients in the Neobianacid group had treatment-related AEs or discontinued the study due to AEs.

In conclusion, starting from the first day of treatment, Neobianacid proved to be noninferior to omeprazole in the relief of heartburn and epigastric pain and burning in the initial 2 weeks, and even better on demand (on average 2–3 tablets/d) compared with omeprazole in the subsequent 2 weeks. In addition, Neobianacid on demand counteracted the predictable worsening of symptoms that follows the suspension of PPI treatment. Furthermore, Neobianacid is a 100% natural product and therefore biodegradable, with no impact on the environment, and it showed high safety without affecting the gut microbiota.

- **Clinical trial: Efficacy and safety of Mipolixin compared to Poliprotect® in the relief of symptoms of moderate functional dyspepsia and heartburn: a randomized, double-blind, parallel group, non-inferiority clinical study.**

This was a multi-centre, randomized, double-blind parallel group, non-inferiority clinical study that evaluated the clinical efficacy of Poliprotect (Neobianacid 1,55g) compared to Mipolixin (Neobianacid's equivalent medical device manufactured by Aboca) in improving overall symptom severity in patients with symptoms of functional dyspepsia and/or heartburn of moderate severity.

Study objectives: Primary objective of this study was to assess whether Neobianacid and Mipolixin have a comparable symptomatic effect on overall symptom severity in patients with symptoms of functional dyspepsia and/or heartburn of moderate severity.

Secondary objectives evaluated the clinical efficacy of both products at different time points during the study, based on improvement of overall symptom severity and specific moderate symptoms of dyspepsia (bothersome postprandial fullness, bothersome early satiation, bothersome epigastric pain, bothersome epigastric burning) and/or bothersome heartburn.

They also evaluated the patient satisfaction and compliance with treatment in both study arms, as well as its safety and tolerability.

Inclusion and exclusion criteria:

- Were included: Adult patients with symptoms of functional dyspepsia of moderate severity and/or bothersome heartburn within 2 and 14 days prior to the screening/baseline visit, with the ability to provide informed consent. According to the current Rome IV diagnostic criteria, functional dyspepsia is

defined as one or more of the following symptoms: postprandial fullness (classified as postprandial distress syndrome), early satiation (inability to finish a normal sized meal, also classified as postprandial distress syndrome), and epigastric pain or burning (classified as epigastric pain syndrome). Patients had to agree not to alter their diet or lifestyle for the duration of the trial. They could also have stable chronic therapies, if not related to the gastroenterological pathologies object of the study, according to the opinion of the investigator.

- Were excluded: Patients presenting at least one of the symptoms or conditions among those listed in the protocol (e.g., anaemia, chronic gastrointestinal bleeding, anorexia, persistent or recurrent vomiting, dysphagia or odynophagia; malignancy and infectious disease; severe heart, pulmonary, liver or kidney disease; gastrointestinal disorders at screening; mental or metabolic disorders; patients receiving any treatment (pharmacological or medical device) for heartburn or dyspeptic symptoms within the last 14 days prior to randomization or any drugs that could affect symptoms or the study; known hypersensitive or intolerance to any components of the study medical devices; pregnant or nursing women).

Study design and methodology: The performance of the study included a total of 2 in-person visits plus 3 scheduled telephone contacts during a period of 2 weeks as follows: a baseline visit (screening and randomization) (day 0), 2 telephone calls at 3 and 7 days after treatment initiation, a final visit at 2 weeks (day 15, 14 days after treatment initiation) and a follow-up telephone call for safety purposes 14 days after the last administration of the study treatment.

A total of 157 patients were 1:1 randomized to receive study treatment: 82 in Neobianacid arm and 75 in Mipolixin arm. The demographic characteristics of patients did not significantly differ between patients receiving Mipolixin and those taking Neobianacid. Overall, 75% and 62.8% of patients treated with Mipolixin and Neobianacid were female and the mean age was 45 years.

The two arms were submitted to 2 different regimens:

- Arm A (experimental arm): Mipolixin one 1.55g chewable tablet 5 times a day for 2 weeks.
- Arm B (active comparator): Neobianacid one 1.55 g chewable tablet 5 times a day for 2 weeks.

The appearance, consistency, packing and labelling of Neobianacid and Mipolixin were identical in order to maintain the double-blind nature of the study.

Endpoints: The primary endpoint was the change in the score of a visual analogue scale (VAS) (from “no symptoms” to “overwhelming symptoms”) used for evaluating overall symptom severity from baseline (day 0) to day 14 after treatment initiation, between the two study arms. A decrease in VAS score of at least 30% was considered as a clinically meaningful improvement, and therefore subjects achieving a decrease $\geq 30\%$ in this scale would be considered as responders.

Among the secondary endpoints, there was the comparison between Neobianacid and Mipolixin in the mean change of the score of a VAS used for evaluating overall symptom severity and individual specific symptoms at different time points during the study.

Patient level of satisfaction with treatment and treatment compliance were also assessed. Incidence and type of adverse events (AEs) and serious adverse events (SAEs) were evaluated, too.

Summary of results: Neobianacid and Mipolixin showed high and similar response rates, of 71.8% and 72.9%, respectively, with response defined as a decrease in the VAS symptoms severity score $\geq 30\%$, thus demonstrating equal efficacy in improving the severity of symptoms of functional dyspepsia of both clinical subtypes (painful and painless) and/or heartburn after 14 days of treatment.

In conclusion, both medical devices demonstrated a high improvement rate of symptoms of over 70%. The efficacy of Neobianacid and Mipolixin has been evident in terms of improvement of each symptom of dyspepsia (postprandial fullness, early satiation, epigastric pain, epigastric burning) and heartburn. VAS score change over time for all those symptoms was found to be statistically significant. No significant differences were observed in the VAS score change over time between the two study arms.

This study has also demonstrated a high rate of treatment compliance and satisfaction with treatment.

None of the patients experienced SAEs or discontinued the study treatment because of AEs. Safety data suggest that Neobianacid is a tolerable therapeutic option in patients with dyspepsia and/or heartburn.

5.3.2 Surveys

In addition to the clinical trials described above, as part of the post-market surveillance activities, different surveys were conducted on Neobianacid, assessing the product in a real-world setting.

Overall, thousands of questionnaires addressed to doctors (both general Practitioners (GPs) and specialists), consumers and pharmacists were assessed, focusing on the perceived effectiveness, tolerability, safety and on the pattern of usage of Neobianacid (tablets, sachets and paediatric).

The surveys aimed at the collection and processing of information to evaluate and/or substantiate the efficacy and safety profile of the product, but also to assess the satisfaction level associated with the use of the product, the quality of life after treatment and to identify possible systematic misuse or off-label use of the device.

Some of these questionnaires available on Neobianacid were performed through an electronic validated platform². That platform is meant to provide with a continuous and proactive collection of data by allowing the users to answer to online surveys based on specific questionnaires. The questionnaires, for each specific category of participants (consumers, physicians, and pharmacists) were developed by clinical experts.

Almost all questions are closed multiple choice questions. The survey participants are also requested to spontaneously report undesirable side effects or product interactions with concomitant treatments addressing any safety issues directly to the Aboca Vigilance Department at a specific mailbox.

Before submitting the questionnaires, a repeatability study of the questionnaires was conducted to assess the degree of precision expressed by the survey tools.

Overall, the surveys for Neobianacid from all doctors and pharmacists consistently reported a prevalent usage of this product to treat the symptoms of gastroesophageal reflux and painful dyspepsia, such as retrosternal pyrosis (heartburn), acid regurgitation, irritative cough and epigastric pain. A large fraction of consumers used the product for digestive difficulties as well (non-painful dyspepsia), such as heaviness, early satiation and bloating.

The high level of agreement across all questionnaires supported the evidence of Neobianacid as a device associated with positive performance, rapidity in symptom resolution/improvement and high level of safety and tolerability, even in special populations such as children over 1 year of age, pregnant and breastfeeding women, with a positive benefit/risk ratio as coherently described by consumers, pharmacists and physicians.

Neobianacid was mainly prescribed as a stand-alone treatment, but there's a consistent number of physicians reporting to have prescribed the product in combination or alternation with anti-secretory drugs.

Furthermore, a high percentage of pharmacists and physicians reported satisfied with the use of Neobianacid in prevention, in situations that cause irritation of the mucosa (e.g., taking anti-inflammatory drugs, pain relievers, periods of mental and physical stress).

Overall, no relevant findings about adverse events and off label use have been detected.

Part of the data of the surveys above mentioned were reported and discussed in a scientific paper published³ in a peer-reviewed journal as a cross-sectional research article on the safety and effectiveness of Neobianacid in the treatment of common functional upper gastrointestinal disorders.

² Cioeta R, Cossu A, Giovagnoni E, Rigoni M and Muti P (2022), A new platform for post-marketing surveillance and real-world evidence data collection for substance-based medical devices. *Front. Drug. Saf. Regul.* 2:992359.

³ Roberto Cioeta, Paola Muti, Marta Rigoni, Luigi Morlando, Filippo Siragusa, Andrea Cossu and Emiliano Giovagnoni. Effectiveness and tolerability of Poliprotect, a natural mucosal protective agent for gastroesophageal reflux disease and dyspepsia: Surveys from patients, physicians, and pharmacists. *Front. Drug Saf. Regul., 2022. Sec. Substance-Based Medical Devices; 2: 1-11*

This study is based on clinical observational data collected from 4,467 subjects (GPs and medical specialists, pharmacists and patients), focusing on the perceived effectiveness, tolerability, safety, and on the pattern of usage of Neobianacid. This study also describes a new methodology of validated surveys that have been developed to fulfil the requirements of the new EU regulation for Substance-based medical devices (SBMD). In conclusion, based on all the surveys available on the product, Neobianacid can be considered a safe and effective product in the treatment of upper gastrointestinal disorders as indicated in the labelling in a large array of subjects, including children and pregnant women.

5.3.3 Other clinical data

The safety data collected from the Vigilance system of Aboca confirm that there have been no serious incidents or safety issues related to the use of the product.

Since no serious incidents have been registered, no cases were reported to Competent Authority because they did not match the criteria of reportability. Furthermore, no trend or safety signal have been identified.

Therefore, it can be concluded that the use of the Neobianacid medical device can be considered safe.

5.4 An overall summary of the clinical performance and safety

Neobianacid is indicated for the treatment of gastroesophageal reflux, stomach pain and digestive difficulties.

The product is also suitable for use during pregnancy, breastfeeding and in children.

The safe and effective profile of the medical device was demonstrated by the post-market clinical follow-up activities and investigations reported above, consisting in two randomized clinical trials and several large-scale surveys. Such evidence is supported also with Vigilance data, as well as pre-clinical testing and assessments.

In all clinical trials the product proved effective and safe in its intended use.

Such evidence from interventional studies is overall coherent to the results obtained with the surveys, which demonstrate a efficacy and safety of the product for all its intended uses in a real-world setting, as directly reported by both healthcare professionals and consumers.

The findings from the surveys also reveal positive results about the efficacy and safety profile of the product in children, in pregnant and breastfeeding women. No safety issues or serious incidents with Neobianacid were reported in this specific populations.

The safety profile of the product is further confirmed by the Vigilance Data recorded for Neobianacid (with referral to more than 25 million packs of both tablets and sachets commercialized worldwide to March 2022), where no serious incidents or safety issues related to the use of the product were identified.

It can be concluded that, based on the analysis of the clinical (and preclinical) data, the benefit/risk ratio of the product is favorable and positive.

5.5 Ongoing or planned post-market clinical follow-up

In the year 2023 other activities have been planned for Neobianacid with the aim to gather information, monitor and follow-up on the clinical efficacy and safety of the product.

Those activities consist on the survey conducted through an electronic validated platform (cited in section 5.3) on the experience of the use of the product, collected by physicians, pharmacists and patients.

The safety profile of the product will be continuously monitored through the collection and analysis of Vigilance data.

6. POSSIBLE DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES

In general, available treatments used in gastroesophageal reflux and dyspeptic symptoms are listed below⁴:

- Proton pump inhibitors (PPIs)
- H2 Blockers
- Antacids
- Alginates
- Sucralfate
- Prokinetics

Different drugs, such as PPIs, are highly effective in preventing complications and reducing symptoms⁵ of acid-induced esophageal and gastroduodenal diseases. Proton pump inhibitors are considered the gold standard treatment for gastroesophageal reflux disease (GERD) and functional dyspepsia⁶.

Neobianacid can represent a valid therapeutic option as it is effective for the treatment of both dyspepsia and gastro-oesophageal reflux disease, acting by enhancing and strengthening the physiological defense systems without intervening directly on acid secretion. This clinical benefit is expressly due to the non-pharmacological mode of action of the whole product⁷, which keeps intact (and enhances) the natural mucosal defense system.

Nonpharmacological gastro-esophageal disorders management is the gap filled by Neobianacid, and the innovation brought forward by the product.

7. SUGGESTED PROFILE AND TRAINING FOR USERS

The device is intended for lay users.

⁴ L. Wauters, R. Dickman, V. Drug, A. Mulak, J. Serra, P. Enck, J. Tack, E.F.C. Group, A. Accarino, G. Barbara, S. Bor, B. Coffin, M. Corsetti, H. De Schepper, D. Dumitrascu, A. Farmer, G. Gourcerol, G. Hauser, T. Hausken, G. Karamanolis, D. Keszthelyi, C. Malagelada, T. Milosavljevic, J. Muris, C. O'Morain, A. Papathanasopoulos, D. Pohl, D. Romyantseva, G. Sarnelli, E. Savarino, J. Schol, A. Sheptulin, A. Smet, A. Stengel, O. Storonova, M. Storr, H. Tornblom, T. Vanuytsel, M. Velosa, M. Waluga, N. Zarate, and F. Zerbib, United European Gastroenterology (UEG) and European Society for Neurogastroenterology and Motility (ESNM) consensus on functional dyspepsia. *United European Gastroenterol J* 9 (2021) 307-331.

⁵ B. Scally, J.R. Emberson, E. Spata, C. Reith, K. Davies, H. Halls, L. Holland, K. Wilson, N. Bhala, C. Hawkey, M. Hochberg, R. Hunt, L. Laine, A. Lanas, C. Patrono, and C. Baigent, Effects of gastroprotectant drugs for the prevention and treatment of peptic ulcer disease and its complications: a meta-analysis of randomised trials. *Lancet Gastroenterol Hepatol* 3 (2018) 231-241

⁶ D. Keszthelyi, C. Malagelada, T. Milosavljevic, J. Muris, C. O'Morain, A. Papathanasopoulos, D. Pohl, D. Romyantseva, G. Sarnelli, E. Savarino, J. Schol, A. Sheptulin, A. Smet, A. Stengel, O. Storonova, M. Storr, H. Tornblom, T. Vanuytsel, M. Velosa, M. Waluga, N. Zarate, and F. Zerbib, United European Gastroenterology (UEG) and European Society for Neurogastroenterology and Motility (ESNM) consensus on functional dyspepsia. *United European Gastroenterol J* 9 (2021) 307-331

⁷ The concept of non-pharmacological mechanism of action in medical devices made of substances in practice: what pharmacology can do to promote the scientific implementation of the European medical device regulation *Pharmadvances* 2020 01s: 4-12 doi:10.36118/pharmadvances.01.2020.02s

8. REFERENCE TO ANY HARMONISED STANDARDS AND CS APPLIED

- ISO 14971: current version “Medical devices — Application of risk management to medical devices.
- ISO 10993-1: current version “Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process”
- ISO 15223-1: current version “Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1 – General Requirements”
- ISO 9001: current version “Quality management system”
- ISO 13485: current version “Medical devices – Quality management systems – Requirements for regulatory purposes”.

9. REVISION HISTORY

| SSCP revision number | Date issued | Change description | Revision validated by the Notified Body |
|-----------------------------|--------------------|---------------------------|--|
| 1 | 24/06/2022 | First issue | |
| 2 | 11/09/2023 | Data update | Yes. Validation language: English |

[Break page]

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A summary of the safety and clinical performance of the device, intended for patients, is given below.

Document revision: 2

Date issued: 11/09/2023

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace the Instructions For Use to provide information on the safe use of the device.

1 DEVICE IDENTIFICATION AND GENERAL INFORMATION

○ Device trade name

NEOBIANACID, NEOBIANACID PEDIATRIC

○ Manufacturer; name and address

Aboca S.p.A. Società Agricola - Loc. Aboca 20 – 52037 Sansepolcro (AR) – Italy

○ Basic UDI-DI

8032472BIOANAbY2AA

○ Year when the device was first CE-marked

2014 (Certification according to Directive 93/42/EEC)

2 INTENDED USE OF THE DEVICE

○ Intended purpose

Neobianacid is indicated for the treatment of gastroesophageal reflux and stomach pain (painful dyspepsia), characterized by acid regurgitation, burning, pain, irritative cough, dysphonia.

Neobianacid is also indicated in the treatment of symptoms associated with digestive difficulties (non-painful dyspepsia), such as acidity, heaviness, and bloating after meals, both recurrent and in cases of copious meals.

Neobianacid can be useful in prevention, in situations that cause irritation of the mucosa (e.g., taking anti-inflammatory drugs, periods of mental and physical stress).

The product is also suitable for use in pregnancy, during lactation and in children over 6 years of age.

Neobianacid Pediatric is indicated from one year of age for the treatment of gastroesophageal reflux and acidity characterized by regurgitation, vomiting, irritative cough and abdominal pain.

Neobianacid Pediatric is also indicated in the treatment of symptoms associated with dyspepsia (digestive difficulty), such as pain and heartburn, heaviness and bloating after meals.

○ Indications and intended patient groups

Neobianacid is indicated in adults and children over 6 years of age. The device is also suitable for use in pregnancy and during breastfeeding.

Neobianacid Pediatric is indicated in children over 1 year of age.

- **Contraindications**

The device shall not be used in case of individual hypersensitivity or allergy to one or more components.

However, if symptoms persist, consult your doctor or pharmacist.

It is recommended to take the device at least two hours before or after taking medicines, particularly tetracyclines. There are no known interactions of Neobianacid or Neobianacid Pediatric with other medical devices or foods.

3 DEVICE DESCRIPTION

- **Device description and material/substances in contact with patient tissues**

Neobianacid is a medical device made of substance and it is available in tablets and granular sachets. Neobianacid Pediatric is available in granular sachets.

Neobianacid needs to be dissolved in the mouth and Neobianacid Pediatric needs to be dissolved in little water and then ingested.

Neobianacid and Neobianacid pediatric are 100% natural, biodegradable and with organic ingredients. They contain Poliprotect, a vegetal and mineral complex. They do not contain synthetic, semi-synthetic or genetically modified substances and the manufacturing processes with which it is obtained do not use solvents or other chemically synthesized substances.

The qualitative composition of Neobianacid is described below:

One 1.55 g tablet or sachet contains:

- Poliprotect® natural complex 528 mg, based on
 - Polysaccharides (LPME extractive fractions of Mallow* and Marshmallow*; dehydrated Aloe vera leaf gel* - 105 mg) and Minerals and
 - Minerals (Limestone and Nahcolite 423 mg),
- LPME extractive flavonoid fraction of Chamomile* and Liquorice*, - 46 mg;

Brown Sugar*, Arabic Gum, Natural mint flavour or Natural lemon flavour.

*Ingredient from organic farming.

Each tablet contains 136 mg (3.4 mmol) of calcium.

The qualitative composition of Neobianacid pediatric is described below:

A 775 mg sachet of Neobianacid pediatric contains:

- Poliprotect® natural complex 264 mg, based on
 - Polysaccharides (LPME extractive fractions of Mallow* and Marshmallow*; dehydrated Aloe vera leaf gel* - 52 mg) and
 - Minerals (Limestone and Nahcolite 212 mg);
- LPME extractive flavonoid fraction of Chamomile* and Liquorice* - 23 mg;

Brown sugar*; Arabic gum; Natural flavour.

*Ingredient from organic farming.

Each sachet contains 68 mg (1.7 mmol) of calcium.

In both Neobianacid and Neobianacid pediatric, the extractive fractions contained in the product are obtained through an exclusive technological process, called LPME (Liquid-Phase Microextraction), developed by Aboca research.

- **Information about medicinal substances in the device, if any**

Not applicable – the device does not contain medicinal substances

- **Description of how the device is achieving its intended mode of action**

NeoBianacid enhances the protective action of the mucus that naturally covers the mucosa and quickly counteracts heartburn, pain and a feeling of acidity.

- **Description of accessories, if any**

Not applicable – no accessories are intended to be used in combination with NeoBianacid

4 RISKS AND WARNINGS

Contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

- **How potential risks have been controlled or managed**

The manufacturer manages the risk analysis and risk control according to the relevant harmonised standards.

- **Remaining risks and undesirable effects**

The main residual risks identified for NeoBianacid are the following:

- Subjective allergic reaction of the user – although the product is non-sensitising, it is possible that a small number of people may be allergic to one or more components of the product; the labelling and the instructions for use list the components of the product, so the user has all the information needed in order to verify the ingredients and therefore the safe use in cases of allergies.
- Microbial contamination in case of poor storage and/or use of the product – in order to minimize this risk, the storage conditions of the product are clearly reported on the labelling and instructions for use.

NeoBianacid is normally well tolerated, however, if the user notices an adverse effect potentially attributable to the product or its use, it is recommended to inform a doctor and/or pharmacist and the Manufacturer Aboca.

- **Warnings and precautions**

Do not use in case of individual hypersensitivity or allergy to one or more components. If symptoms persist, seek medical advice. The glycyrrhizin content from the Liquorice in NeoBianacid is lower than that reported to cause a possible increase in blood pressure.

- **Summary of any field safety corrective action, (FSCA including FSN) if applicable**

No FSCA has been implemented for the device NeoBianacid since its first date of commercialization in 2014.

5 SUMMARY OF CLINICAL EVALUATION AND POST-MARKET CLINICAL FOLLOW-UP

- **Background of the device**

NeoBianacid is a 100% natural, gluten-free, lactose-free medical device indicated for the treatment of gastroesophageal reflux, painful dyspepsia and digestive difficulties.

The product is also suitable for use during pregnancy, breastfeeding and in children.

NeoBianacid and NeoBianacid pediatric_act rapidly forming a barrier effect film with the following activities:

PROTECTIVE - thanks to the high affinity for the gastroesophageal mucosa, it forms a film that protects and promotes the physiological reparative processes.

LOCAL NEUTRALIZING OF ACIDITY - strengthens the local neutralizing properties of gastric mucus, hindering the contact of acids with the gastric mucosa and prevents damage

ANTIOXIDANT - counteracts irritation due to free radicals.

This way Neobianacid and Neobianacid pediatric protect the mucosa from the contact with gastric juices and irritating substances.

Neobianacid and Neobianacid pediatric modulate the acidity without suppressing it, thus respecting the physiology of the stomach and the functionality of the intestinal flora.

- **The clinical evidence of Neobianacid on safety and efficacy**

The efficacy and safety profile of Neobianacid is demonstrated by clinical data coming from:

- two randomized clinical trials
- several large-scale surveys in a real-world setting addressed to both healthcare professionals and consumers
- safety data collected from the manufacturer’s Vigilance system

About the clinical trials, one of them was a non-inferiority trial ¹ conducted on Neobianacid to demonstrate its clinical efficacy and safety compared to Omeprazole 20 mg/day in the treatment of symptoms of non-erosive gastroesophageal reflux disease and painful dyspepsia.

The study demonstrated that the efficacy of Neobianacid (Poliprotect), assessed on days 14 (V1), 28 (V2) and 56 (V3) from the start of treatment, is comparable to that of omeprazole 20 mg (PPI) in the reduction of symptoms associated with gastroesophageal reflux and painful dyspepsia, as shown in Figure 1. This efficacy was assessed using a validated measurement tool (VAS, rating scale from 0 to 100, where 0 indicated "no symptoms" and 100 indicated "overwhelming symptoms"), by which each patient reported the intensity of the symptoms.


Figure 1.

Neobianacid was well tolerated and showed an excellent benefit / risk ratio, improving overall symptoms and quality of life. Furthermore, it showed no effect on microbiota composition.

In conclusion, starting from the first day of treatment, Neobianacid proved noninferior to omeprazole in the relief of heartburn and epigastric pain and burning in the initial 2 weeks, and even better on demand compared with omeprazole in the subsequent 2 weeks. In addition, Neobianacid on demand counteracted the predictable worsening of symptoms that follows the suspension of PPI treatment.

¹ Enrico Stefano Corazziari et al. Poliprotect vs Omeprazole in the Relief of Heartburn, Epigastric Pain, and Burning in Patients Without Erosive Esophagitis and Gastroduodenal Lesions: A Randomized, Controlled Trial. *Am J Gastroenterol* 2023;00:1–11. <https://doi.org/10.14309/ajg.0000000000002360>

Unlike with Neobianacid, the initial benefit of omeprazole is maintained against significantly higher usage of rescue medicine compared with the Neobianacid group and despite daily omeprazole intake.

The other clinical trial aimed to evaluate the clinical efficacy of Neobianacid compared to its equivalent medical device Mipolixin (manufactured by Aboca) in improving the overall symptom severity of functional dyspepsia and/or heartburn.

From the analysis of results, both products enable a high improvement rate of overall severity of symptoms of functional dyspepsia of both clinical subtypes (painful and painless) and/or heartburn, showing a comparable efficacy. This study has also demonstrated a high rate of treatment compliance and satisfaction with treatment, as well as confirming the safety profile and excellent tolerability of the product.

In addition to the clinical trials described above, the surveys conducted on Neobianacid (available as tablets, sachets and paediatric formulation) represent a high level of demonstration of the clinical and safety profile of the product in a real-world setting, directly gathered from both healthcare professionals and consumers.

The surveys aimed at the collection of information to evaluate the efficacy and safety profile of the product, but also to assess the satisfaction level associated with the use of the product, the quality of life after treatment and to identify possible systematic misuse of the device.

Some surveys were performed through an electronic validated platform². That platform is meant to provide with a continuous collection of data by allowing the users to answer to online surveys based on specific questionnaires.

Overall, the thousands of surveys from doctors, herein included different specialists, pharmacists and consumers consistently reported a prevalent usage of this product to treat the symptoms of gastroesophageal reflux and painful dyspepsia, such as retrosternal burning (heartburn), acid regurgitation, stomach pain, irritative cough. A large fraction of consumers used the product for digestive difficulties as well, such as heaviness, early satiation and bloating.

The high level of agreement across all questionnaires supports the evidence of Neobianacid as a device associated with positive assessment of its effectiveness, high speed of symptom resolution/improvement and safety and tolerability with a positive benefit/risk ratio as concordantly described by consumers, pharmacists and physicians.

Part of the data of the surveys above mentioned were reported and discussed in a scientific paper published³ in a peer-reviewed journal as a cross-sectional research article on the safety and effectiveness of Neobianacid in the treatment of common functional upper gastrointestinal disorders.

The surveys confirmed the safe and effective profile of the product in the paediatric population and even during pregnancy and breastfeeding.

The safety profile of the product in all its intended users is further confirmed by the Vigilance Data recorded for Neobianacid through the Vigilance system **Errore. Il segnalibro non è definito.** (with referral to more than 25 million packs of both tablets and sachets commercialized worldwide to March 2022), where no serious incidents or safety issues related to the use of the product were identified.

In the year 2023 other activities have been planned for Neobianacid with the aim to gather information, monitor and follow-up on the clinical efficacy and safety of the product.

Furthermore, the safety profile of the product will be continuously monitored through the collection and analysis of Vigilance data.

² Cioeta R, Cossu A, Giovagnoni E, Rigoni M and Muti P (2022), A new platform for post-marketing surveillance and real-world evidence data collection for substance-based medical devices. *Front. Drug. Saf. Regul.* 2:992359

³ Roberto Cioeta, Paola Muti, Marta Rigoni, Luigi Morlando, Filippo Siragusa, Andrea Cossu and Emiliano Giovagnoni. Effectiveness and tolerability of Poliprotect, a natural mucosal protective agent for gastroesophageal reflux disease and dyspepsia: Surveys from patients, physicians, and pharmacists. *Front. Drug Saf. Regul., 2022. Sec. Substance-Based Medical Devices; 2: 1-11*

6 POSSIBLE DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES

When considering alternative treatments, it is recommended to contact your healthcare professional who can take into account your individual situation.

○ General description of therapeutic alternatives

In general, available treatments used in gastroesophageal reflux and dyspeptic symptoms are listed below ⁴:

- Proton pump inhibitors (PPIs)
- H2 Blockers
- Antacids
- Alginates
- Sucralfate
- Prokinetics

Different drugs, such as PPIs, are highly effective in preventing complications and reducing symptoms⁵ of acid-induced esophageal and gastroduodenal diseases. Proton pump inhibitors are considered the gold standard treatment for gastroesophageal reflux disease (GERD) and functional dyspepsia ⁶.

Neobianacid can represent a valid therapeutic option as it is effective for the treatment of both dyspepsia and gastro-oesophageal reflux disease, acting by enhancing and strengthening the physiological defense systems without intervening directly on acid secretion. This clinical benefit is expressly due to the non-pharmacological mode of action of the whole product⁷, which keeps intact (and enhances) the natural mucosal defense system.

Nonpharmacological gastro-esophageal disorders management is the gap filled by Neobianacid, and the innovation brought forward by the product.

7 SUGGESTED TRAINING FOR USERS

The device is intended for lay users. No training is required.

⁴ L. Wauters, R. Dickman, V. Drug, A. Mulak, J. Serra, P. Enck, J. Tack, E.F.C. Group, A. Accarino, G. Barbara, S. Bor, B. Coffin, M. Corsetti, H. De Schepper, D. Dumitrascu, A. Farmer, G. Gourcerol, G. Hauser, T. Hausken, G. Karamanolis, D. Keszthelyi, C. Malagelada, T. Milosavljevic, J. Muris, C. O'Morain, A. Papathanasopoulos, D. Pohl, D. Romyantseva, G. Sarnelli, E. Savarino, J. Schol, A. Sheptulin, A. Smet, A. Stengel, O. Storonova, M. Storr, H. Tornblom, T. Vanuytsel, M. Velosa, M. Waluga, N. Zarate, and F. Zerbib, United European Gastroenterology (UEG) and European Society for Neurogastroenterology and Motility (ESNM) consensus on functional dyspepsia. *United European Gastroenterol J* 9 (2021) 307-331.

⁵ B. Scally, J.R. Emberson, E. Spata, C. Reith, K. Davies, H. Halls, L. Holland, K. Wilson, N. Bhala, C. Hawkey, M. Hochberg, R. Hunt, L. Laine, A. Lanas, C. Patrono, and C. Baigent, Effects of gastroprotectant drugs for the prevention and treatment of peptic ulcer disease and its complications: a meta-analysis of randomised trials. *Lancet Gastroenterol Hepatol* 3(2018) 231-241.

⁶ D. Keszthelyi, C. Malagelada, T. Milosavljevic, J. Muris, C. O'Morain, A. Papathanasopoulos, D. Pohl, D. Romyantseva, G. Sarnelli, E. Savarino, J. Schol, A. Sheptulin, A. Smet, A. Stengel, O. Storonova, M. Storr, H. Tornblom, T. Vanuytsel, M. Velosa, M. Waluga, N. Zarate, and F. Zerbib, United European Gastroenterology (UEG) and European Society for Neurogastroenterology and Motility (ESNM) consensus on functional dyspepsia. *United European Gastroenterol J* 9 (2021) 307-331

⁷ The concept of non-pharmacological mechanism of action in medical devices made of substances in practice: what pharmacology can do to promote the scientific implementation of the European medical device regulation *Pharmadvances* 2020 01s: 4-12 doi:10.36118/pharmadvances.01.2020.02s

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1.3 Tootja single registration number (SRN)

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⁵ B. Scally, J.R. Emberson, E. Spata, C. Reith, K. Davies, H. Halls, L. Holland, K. Wilson, N. Bhala, C. Hawkey, M. Hochberg, R. Hunt, L. Laine, A. Lanas, C. Patrono, and C. Baigent, Effects of gastroprotectant drugs for the prevention and treatment of peptic ulcer disease and its complications: a meta-analysis of randomised trials. Lancet Gastroenterol Hepatol 3 (2018) 231-241

⁶ D. Keszthelyi, C. Malagelada, T. Milosavljevic, J. Muris, C. O'Morain, A. Papathanasopoulos, D. Pohl, D. Romyantseva, G. Sarnelli, E. Savarino, J. Schol, A. Sheptulin, A. Smet, A. Stengel, O. Storonova, M. Storr, H. Tornblom, T. Vanuytsel, M. Velosa, M. Waluga, N. Zarate, and F. Zerbib, United European Gastroenterology (UEG) and European Society for Neurogastroenterology and Motility (ESNM) consensus on functional dyspepsia. United European Gastroenterol J 9 (2021) 307-331

⁷ The concept of non-pharmacological mechanism of action in medical devices made of substances in practice: what pharmacology can do to promote the scientific implementation of the European medical device regulation Pharmadvances 2020 01s: 4-12 doi:10.36118/pharmadvances.01.2020.02s

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